



# American Access

## Portable Ramp Selection Guidelines

To select your ramp, first determine the rise for your application.  
Never exceed a slope greater than 2" on 12" with an occupied chair or scooter.  
Never exceed a slope greater than 3" on 12" with an unoccupied chair or scooter

Ramp Length	ADA recommendations (1:12 Slope) Unassisted	Maximum Rise Occupied (2:12 Slope) Assisted	Maximum Rise Unoccupied (3:12 Slope) Unoccupied
1 foot	1"	2"	3"
2 feet	2"	4"	6"
3 feet	3"	6"	9"
4 feet	4"	8"	12"
5 feet	5"	10"	15"
6 feet	6"	12"	18"
7 feet	7"	14"	21"
8 feet	8"	16"	24"
9 feet	9"	18"	27"
10 feet	10"	20"	30"
11 feet	11"	22"	33"
12 feet	12"	24"	36"

American Access, Inc. \* 8380 Wolf Lake Drive, Suite 108 \* Bartlett, TN 38133 \* Phone 888-790-9269  
www.americanaccess.net \* email: sales@americanaccess.net